

## Data Sheet

Page 1/1



### 30M X 0.32MM ID - BPX5 0.25UM GC Column

**Part Number:** TJ-054145  
**Manufacturer:** Trajan Scientific  
**Quantity:** 1 pc

## Description

30M X 0.32MM ID - BPX5 0.25UM GC-Säule

## Technical Information

|                                  |              |                                  |                                         |
|----------------------------------|--------------|----------------------------------|-----------------------------------------|
| <b>Brand:</b>                    | BP           | <b>Base Material:</b>            | 5% Phenyl 95% polysilphenylene siloxane |
| <b>Product Type:</b>             | GC Column    | <b>Temperature Minimum [°C]:</b> | -40                                     |
| <b>Length:</b>                   | 30 m         | <b>Temperature Maximum [°C]:</b> | 360                                     |
| <b>Internal Diameter:</b>        | 0.32 mm      | <b>Mode of Separation:</b>       | Gas Chromatography                      |
| <b>Film Thickness [µm]:</b>      | 0.25 µm      |                                  |                                         |
| <b>Chemistry:</b>                | 5 (5HT, 5MS) |                                  |                                         |
| <b>Chemistry Producers Term:</b> | BPX5         |                                  |                                         |
| <b>USP Class:</b>                | G27, G36     |                                  |                                         |
| <b>Country of Origin:</b>        | Australia    |                                  |                                         |