

Data Sheet

Page 1/1



60M X 0.32MM ID - BP5 0.25UM GC Column

Part Number: TJ-054178
Manufacturer: Trajan Scientific
Quantity: 1 pc

Description

60M X 0.32MM ID - BP5 0.25UM GC-Säule

Technical Information

Brand:	BP	Base Material:	5% Phenyl 95% methylpolysiloxane
Product Type:	GC Column	Temperature Minimum [°C]:	-60
Length:	60 m	Temperature Maximum [°C]:	320
Internal Diameter:	0.32 mm	Mode of Separation:	Gas Chromatography
Film Thickness [µm]:	0.25 µm		
Chemistry:	5 (5HT, 5MS)		
Chemistry Producers Term:	BP5		
USP Class:	G27, G36		
Country of Origin:	Australia		