

Data Sheet

Page 1/1



60M X 0.53MM ID - BP20 1.0UM GC Column

Part Number: TJ-0544515
Manufacturer: Trajan Scientific
Quantity: 1 pc

Description

60M X 0.53MM ID - BP20 1.0UM GC-Säule

Technical Information

Brand:	BP	Base Material:	Polyethylene glycol
Product Type:	GC Column	Temperature Minimum [°C]:	-40
Length:	60 m	Temperature Maximum [°C]:	330
Internal Diameter:	0.53 mm	Mode of Separation:	Gas Chromatography
Film Thickness [µm]:	1 µm		
Chemistry:	20		
Chemistry Producers Term:	BP20		
USP Class:	G16, G20		
Country of Origin:	Australia		