

Data Sheet

Page 1/1



60M X 0.32MM SOLGEL-WAX 0.25UM GC Column

Part Number: TJ-054789
Manufacturer: Trajan Scientific
Quantity: 1 pc

Description

60M X 0.32MM SOLGEL-WAX 0.25UM GC-Säule

Technical Information

Brand:	SolGel	Base Material:	Polyethylene glycol
Product Type:	GC Column	Temperature Minimum [°C]:	30
Length:	60 m	Temperature Maximum [°C]:	260
Internal Diameter:	0.32 mm	Mode of Separation:	Gas Chromatography
Film Thickness [µm]:	0.25 µm		
Chemistry:	WAX (WAX MS)		
Chemistry Producers Term:	SOLGEL-WAX		
USP Class:	G16, G35		
Country of Origin:	Australia		