

TQC - Thermimport Quality Control, hereby declares that the instrument mentioned in this declaration has been produced complying with our internal standards and if applicable with the relevant international standards.

The instrument has been tested according the appropriate quality instruction, which is part of TQC's ISO 9001:2008 quality system, that is regularly audited by DNV GL - Business Assurance as the independent national accredited body, which is traceable through certificate number 181638-2015-AQ-NLD-RvA.

CALIBRATION CERTIFICATE No.: VF2100-961	
Instrument Information	Manufacturer: TQC Name: Pyknometer 50ml RVS Serial number: VF2100-961 Complies with: ISO 2811, DIN 53217, ASTM D1475
Calibration Equipment	Name: Ohaus Classic Serial number: 1325a3 Traceable through: mass piece 100g / 8621487A2011
Environmental Conditions	Temperature: 20,5 °C Humidity: 80 % rH Pressure: 1016 mbar
Traceable Information	The accuracy of the calibration equipment used is traceable through National and / or International Standards.

TEST RESULTS				
Equipment	Quantity	Unit	Nominal value	Actual value
1325a3	Mass Empty	g		145,04
1325a3	Mass Full	g		194,88
	Mass Water	g		49,84
	Vol. at Cal. Temp.	cm ³		50,03
	Temp. Cal.	°C		24
	Vol. at Ref. Temp.	cm ³	50	50,02
	Temp. Ref.	°C		20

Issued by: M. J. N. Nadorp


22-08-2017
INSPECTED
THERMIMPORT

TQC B.V.

Molenbaan 19
2908 LL Capelle aan den IJssel
The Netherlands

+31(0) 10 - 79 00 100
+31(0) 10 - 79 00 129

@ info@tqc.eu
www.tqc.eu

B.T.W. (V.A.T.) NL8078.01.744.B01
Bank ABN-AMRO
BIC ABNANL2A

IBAN (EUR) NL91 ABNA 0598 0945 20
IBAN (USD) NL16 ABNA 0598 0973 25
IBAN (GBP) NL10 ABNA 0598 1025 74

Chamber of Commerce / Kamer van Koophandel nr. 29042770.

We supply under sales conditions registered at the Chamber of Commerce in Rotterdam. Complaints within 8 days after delivery.

Leveringsvoorwaarden zoals gedeponeerd bij K.v.K. te Rotterdam. Reclamaties binnen 8 dagen na levering.

We have
switched banks