

TERUMO® Conventional Needles & Syringes



Superior sharpness and smoothness time after time

- Ultra Sharp Needles
- Clear Barrels
- Bold, Precise Scale Markings
- Smooth Plunger Action
- Tamper Evident Peel Packs

When comfort, performance and cost effectiveness are important, choose Terumo. Our ultra-thin wall needles are precisely honed and beveled for ease of penetration and less patient discomfort. Clear syringe barrels with bold scale markings permit quick visual inspection for precise dosage. An integrated plunger stop helps prevent the plunger from being accidentally pulled out. And our tamper evident blister packs result in minimal waste disposal.


Terumo® Standard Hypodermic Needles

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ NN1825R	18G x 1" TW
■ NN1838R	18G x 1½" TW
■ NN1925R	19G x 1" TW
■ NN1938R	19G x 1½" TW
■ NN2025R	20G x 1" UTW
■ NN2038R	20G x 1½" UTW
■ NN2125R	21G x 1" UTW
■ NN2138R	21G x 1½" UTW
■ NN2225R	22G x 1" UTW
■ NN2238R	22G x 1½" UTW
■ NN2325R	23G x 1" UTW
■ NN2338R	23G x 1½" UTW
■ NN2516R	25G x 5/8" UTW
■ NN2525R	25G x 1" UTW
■ NN2538R	25G x 1½" RW
■ NN2613R	26G x 1/2" RW
■ NN2713R	27G x 1/2" RW
■ NN2732R	27G x 1¼" RW
■ NN3013R	30G x 1/2" RW
■ NN3025R	30G x 1" RW

RW - Regular Wall
TW - Thin Wall
UTW - Ultra Thin Wall

All Competitors



- * In Independent laboratory testing, Terumo needles were shown, on average, to be the sharpest among brands tested:
10% sharper than BD (22G - 30G)
16% sharper than Covidien (22G - 27G)

Call your Terumo representative at:

1.800.888.3786

or visit www.TerumoTMP.com

and click on Find Your Local Terumo Rep.

When you think hypodermics...
Think sharp. Think Terumo.



* Data on file. Terumo Medical Products, April 2016.

These products have no components made of natural rubber latex.

Terumo® Conventional Needles & Syringes

1cc Tuberculin Syringes

Graduations: 0.01cc (1cc total)

TERUMO TB Syringes have clear, bold scale markings which assure dosage control and accuracy. The package is tamper evident, for greater protection of the sterile contents.

Packed: 100/box, 10 boxes/case (1,000)

Product Code Description

SS-01T	1cc TB Syringe without needle, with luer slip
SS-01T2516	1cc TB Syringe with 25G x 5/8" removable needle
SS-01T2609	1cc TB Syringe with 26G x 3/8" removable needle
SS-01T2713	1cc TB Syringe with 27G x 1/2" removable needle



Catheter Tip



Luer Lock Tip



Eccentric Luer Tip



Luer Slip Tip

3cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
SS-03L2025	3cc, 20G x 1" UTW
SS-03L2038	3cc, 20G x 1 1/2" UTW
SS-03L2125	3cc, 21G x 1" UTW
SS-03L2138	3cc, 21G x 1 1/2" UTW
SS-03L2225	3cc, 22G x 1" UTW
SS-03L2238	3cc, 22G x 1 1/2" UTW
SS-03L2325	3cc, 23G x 1" UTW
SS-03L2338	3cc, 23G x 1 1/2" UTW
SS-03L2516	3cc, 25G x 5/8" UTW
SS-03L2525	3cc, 25G x 1" UTW

5cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.2cc (6cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code	Description
SS-05L2025	5cc, 20G x 1" UTW
SS-05L2038	5cc, 20G x 1 1/2" UTW
SS-05L2125	5cc, 21G x 1" UTW
SS-05L2138	5cc, 21G x 1 1/2" UTW
SS-05L2225	5cc, 22G x 1" UTW
SS-05L2238	5cc, 22G x 1 1/2" UTW

10cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.2cc (12cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code	Description
SS-10L2025	10cc, 20G x 1" UTW
SS-10L2038	10cc, 20G x 1 1/2" UTW
SS-10L2125	10cc, 21G x 1" UTW
SS-10L2138	10cc, 21G x 1 1/2" UTW

UTW - Ultra Thin Wall

3cc Syringes without Needle

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
SS-03L	3cc Luer Lock Tip
SS-03S	3cc Luer Slip Tip

5cc Syringes without Needle

Graduations: 0.2cc (6cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code	Description
SS-05L	5cc Luer Lock Tip
SS-05S	5cc Luer Slip Tip

10cc Syringes without Needle

Graduations: 0.2cc (12cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code	Description
SS-10L	10cc Luer Lock Tip
SS-10S	10cc Slip Tip
SS-10ES	10cc Eccentric Luer Tip

20cc Syringes without Needle

Graduations: 1.0cc (25cc total)

Packed: 50/box, 4 boxes/case (200)

Product Code	Description
SS-20L2	20cc Luer Lock Tip
SS-20S2	20cc Luer Slip Tip
SS-20ES2	20cc Eccentric Luer Tip

30cc Syringes without Needle

Graduations: 1.0cc (35cc total)

Packed: 25/box, 8 boxes/case (200)

Product Code	Description
SS-30L	30cc Luer Lock Tip
SS-30S	30cc Luer Slip Tip
SS-30ES	30cc Eccentric Luer Tip

60cc Syringes without Needle

Graduations: 1.0cc (60cc total)

Packed: 25/box, 4 boxes/case (100)

Product Code	Description
SS-60L	60cc Luer Lock Tip
SS-60ES	60cc Eccentric Luer Tip
SS-60C	60cc 2oz. Catheter Tip

TerumoTMP.com

These products have no components made of natural rubber latex.

TERUMO and TERUMO are registered trademarks owned by Terumo Corporation, Tokyo, Japan. SurSaver is a registered trademark owned by Terumo Medical Corporation, Elkton, MD, U.S.A. PM-00687



TERUMO MEDICAL PRODUCTS