



SAFETY MEETS QUALITY

Product Specification

FORMAT

Disposable Nitrile Examination Gloves

- Powder free
- White
- Latex free
- Chemical protection EN 374
- Non-sterile



Material

100% Synthetic Nitrile Latex

Performance requirements

Characteristics	Inspection Level	Acceptable Quality Level
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	1.5
Minor Visual defects	G-1	4.0

Dimensions

Number	1001	1002	1003	1004	1005
Size	XS = 5–6	S = 6–7	M = 7–8	L = 8–9	XL = 9–10
Width	80 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	115 ± 10 mm
Length	min. 240 mm	min. 240 mm	min. 240 mm	min. 240 mm	min. 240 mm

Thickness

Measuring Points	Single Wall (mm)
Cuff (distance 25 mm to the beaded edge)	0.08 – 0.13
Palm (in the middle of the palm)	0.10 – 0.15
Fingertip (distance 15 mm to top of middle finger)	0.10 – 0.17

Physical Properties

Standard	Before Aging	After Aging
Tensile strength (N)	min. 7.5	min. 5.5
Elongation at Break (%)	min. 500	min. 400
Modulus of elasticity (MPa)	min. 21.0	min. 16.0

Powder

No powder lubricant added

Protein content

No proteins contained

Cuff

Beaded cuff

Design

Ambidextrous

Surface

Micro-roughened palm

Color

White

Quality assurance

EG construction sample test BGIA – CE 0121

Performance standards

PSA – EN 374 CAT III, EN 388, EN 420
MPG – EN 455 part 1–4
EG-Standard 89/686/EWG
EG-Standard 93/42/EWG, Appendix II.3 + V.3
EN ISO 9001:2008
EN ISO 13485:2003 + AC:2007

Packing Style

100 pcs per dispenser, 10 dispensers per carton

Hint

The above information apply to the standard version.
Subject to change.