

LABORHITS

VACUU-VIEW and VACUU-VIEW extended

Vacuum measuring instruments

Rough vacuum gauge **VACUU-VIEW** - Measuring with highest precision from atmosphere down to 0.1 mbar

Compact design with integrated sensors

Compatible with gauge DCP 3000 and controller CVC 3000

Rough and fine vacuum gauge **VACUU-VIEW extended** with an extended measurement range from atmosphere down to 10^{-3} mbar

Easy to read, illuminated display

Easy, menu driven handling

Chemically resistant and robust



Specifications	VACUU-VIEW	VACUU-VIEW extended
Measuring range	1100 - 0.1 mbar (hPa)	1100 - 0.001 mbar (hPa)
Measuring range	825 - 0.075 torr	825 - 0.001 torr
Measuring principle	ceramic diaphragm sensor with capacitive readout	ceramic diaphragm sensor + ceramic jacketed Pirani sensor
Accuracy of measurement	±1 mbar (hPa)	±15 % of indicated value in the range from 10 to 0.01 mbar (hPa) ±3 mbar for >10 mbar (hPa)
Vacuum connection	KF DN 16 / hose nozzle DN 6/10 mm	
Dimensions (L x W x H)	103 x 62 x 50 mm	
Weight	190 g	
Power supply	100 - 230 V, 50/60 Hz (CEE/CH/UK/AUS/CN)	
Cat. No.	6.268 865	6.268 864