

## CN-15 CABINET

The CN-15 darkroom provides a large effective capacity and UV power intensity unequalled in this field. The darkrooms offer any combination of UV sources, simultaneously or not. Its key features are:

- Extra large capacity
- Black rubber curtain for easy access into the darkroom
- White-light bulb for normal observation
- UV absorber shield to protect the user from UV light
- Removable bottom panel for use with a Vilber Lourmat ETX fluorescent table



CN-15.LC

## CN-6 CABINET

The CN-6 darkroom holds one or two hand-held UV lamps... (VL-6 model) in any of the three following wavelengths: 254, 365 or 312 nm. The darkroom is supplied without lamps and allows different lighting possibilities according to the user's choice. Its key features are:

- Large capacity
- Black rubber curtain for easy access into the darkroom
- UV absorber shield to protect the user from the UV light
- Removable lamps that can be used for hand-held applications



CN-6

## SPECIFICATIONS

Models	Tubes (Watt)	Wavelength (nm)	Intensity at bottom ( $\mu\text{W}/\text{cm}^2$ )	Size - W x D x H (mm)
CN-15.LL	4 x 15-W	365	2 000	505 x 415 x 280
CN-15.CC	4 x 15-W	254	1 750	505 x 415 x 280
CN-15.MM	4 x 15-W	312	2 500	505 x 415 x 280
CN-15.LC	2 x 15-W	365	1 050	505 x 415 x 280
	2 x 15-W	254	900	
CN-15.LM	2 x 15-W	365	1 050	505 x 415 x 280
	2 x 15-W	312	1 300	
CN-15.MC	2 x 15-W	312	1 300	505 x 415 x 280
	2 x 15-W	254	900	
CN-6	Lamps not included			300 x 280 x 240