


Chromatographic columns with frit and NS socket with PTFE Stopcock



Chromatographic columns with frit and NS socket with PTFE Stopcock

Chromatographic columns with sintered frit (porosity: 0), with PTFE stopcock or needle-valve stopcock (bore: 0-2,5 mm), made of DURAN® tubing

Effectiv length (mm)	Dia. (mm)	Capacity (ml)	Socket (NS)	Stopcock (Description)	Pack (qty.)	Cat. No.
100	10	8	14/23	PTFE stopcock	1	5 4507 00
200	10	15	14/23	PTFE stopcock	1	5 4507 01
300	10	23	14/23	PTFE stopcock	1	5 4507 02
200	15	35	14/23	PTFE stopcock	1	5 4507 03
400	20	125	29/32	PTFE stopcock	1	5 4507 04
600	30	430	29/32	PTFE stopcock	1	5 4507 06