

KERN CM 150-1N

Pocket balance with integrated pocket calculator

KERN®

Article codes, product groups & Status information

Model article-code	CM 150-1N
Model-Serie	CM
Assortment	KERN
Product Type	Laboratory balances
Product group	Pocket balances
Product usage type	Main article
EAN code (Model)	4045761027282
Planned start date	2004-1HJ
Customs tariff number	84238129
CE Label	Yes
Country of Origin	CN

Technical data - Primary

Weighing capacity	150 g
Readability	0,1 g

Technical data - Construction

Casing material	Plastic
Dimension housing (W×D×H)	85×130×25 mm
Material weighing plate	Plastic
Weighing surface (d)	70×80 mm

Technical data - Display

Display height	12 mm
Backlight	Yes

Technical data - Power supply

Battery	2×15 V AAA
Operating time	33 h

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 0,2 g
Weighing capacity	150 g
Readability	0,1 g
Recommended adjusting weight	100 g (M3)
Repeatability	0,1 g
Stabilisation time under laboratory conditions	2 s
Tare range	150 g
Warm up time	10 min
Weighing system	Strain gauge
Weighing units	G;tl (Tw);tol

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	205×136×32 mm
Gross weight	0,35 kg
Net weight	0,20 kg

KERN CM 150-1N

Pocket balance with integrated pocket calculator



Shipping method

Parcel service

Packaging component - by weight - cardboard

454 g

Packaging component - by weight - styrofoam

34 g

Delivery time

1 d

Optional Services

DAkKS Calibration

963-127

Functions

Standard



Option

DAkKS