


Conical joint clips type “Safety Clip”



Conical joint clips type “Safety Clip”

For conical ground joints, material: POM, height adjustable by variable gear rods, also for use with sleeves with grip collar, outstanding chemical resistance, operating temperature up to approx. 90 °C

For ground joints (NS)	Colour	Pack(qyt.)	Cat. No.
12/21	blue	10	1 1038 12
14/23	blue	10	1 1038 14
19/26	blue	10	1 1038 19
24/29	blue	10	1 1038 24
29/32	blue	10	1 1038 29
34/35	blue	10	1 1038 34
45/40	blue	10	1 1038 45