

## Dropping funnels, conical, with PTFE plug



Dropping funnels, conical, with PTFE plug

Conical shape with NS socket and NS cone, with PTFE plug, with screwthread retaining device, with inscription label, without scale, material: borosilicate glass 3.3

| Capacity (ml) | Socket (NS) | Cone (NS) | Bore (mm) | Pack (qty.) | Cat. No.  |
|---------------|-------------|-----------|-----------|-------------|-----------|
| 50            | 19/26       | 19/26     | 2,5       | 1           | 4 4728 02 |
| 100           | 19/26       | 19/26     | 2,5       | 1           | 4 4737 02 |
| 250           | 29/32       | 29/32     | 4,0       | 1           | 4 4749 04 |
| 500           | 29/32       | 29/32     | 4,0       | 1           | 4 4758 04 |
| 1000          | 29/32       | 29/32     | 6,0       | 1           | 4 4770 04 |
| 2000          | 29/32       | 29/32     | 6,0       | 1           | 4 4773 04 |