



IMHOFF CONES - "glass" & "plastic"

**"glass" - "W/O stopcock"**

- Designed for sedimentation tests.
- Manufactured from borosilicate glass 3.3 according to DIN 12672.
- With graduations to 1000 ml. Without stopcock.
- Volume graduations and inscriptions are in white enamel.

catalogue number	graduation range	subdivision	pack quantity
043.20.001	0 - 2 ml	0,1 ml	1 piece
	> 2 - 10 ml	0,5 ml	
	> 10 - 40 ml	1,0 ml	
	> 40 - 100 ml	2,0 ml	
	> 100 - 1000 ml	50,0 ml	

**"glass" - "with stopcock"**

- Designed for sedimentation tests.
- Manufactured from borosilicate glass 3.3 according to DIN 12672.
- With graduations to 1000 ml. With stopcock.
- Volume graduations and inscriptions are in white enamel.

catalogue number	graduation range	subdivision	pack quantity
043.20.002	0 - 2 ml	0,1 ml	1 piece
	> 2 - 10 ml	0,5 ml	
	> 10 - 40 ml	1,0 ml	
	> 40 - 100 ml	2,0 ml	
	> 100 - 1000 ml	50,0 ml	

**"plastic" - "SAN"**

- Designed for sedimentation tests.
- Manufactured from crystal clear special polymer SAN.
- With graduations to 1000 ml. Without stopcock.
- Volume graduations and inscriptions are in embossed form.

catalogue number	graduation range	subdivision	pack quantity
043.21.001	0 - 2 ml	0,1 ml	1 piece
	> 2 - 10 ml	0,5 ml	1 piece
	> 10 - 40 ml	1,0 ml	1 piece
	> 40 - 100 ml	2,0 ml	1 piece
	> 100 - 1000 ml	50,0 ml	1 piece

**IMHOFF CONE RACK**

- Manufactured from PMMA/ stainless steel for holding 2 imhoff sedimentation cones.
- They can carry any type of sedimentation cone as glass or plastic, with or without stopcock.
- Compact design allows easy carry of sedimentation cones even when cones are filled.

catalogue number	pack quantity
043.22.001	1 piece

