

KERN ORA 3AB

Refractive index measurement for laboratories and the industry



Article codes, product groups & Status information

Model article-code	ORA 3AB
Model-Serie	ORA-AL
Assortment	Optics
Product Type	Refractometers
Product group	Analogue refractometers
Product usage type	Main article
EAN code (Model)	4045761261051
Successor Model	ORA 3AA
Planned start date	2022-1HJ
Customs tariff number	90275000
Model Status	Yes
Country of Origin	CN

Technical data - Primary

Scales	Brix original wort (beer)
Scale Range	0 % - 32 % 1,000 SG Wort - 1,130 SG Wort

Technical data - Construction

Type	Handheld refractometer
Display scale	Analogue
Dimensions (W×D×H)	198×40×40 mm

Technical data - Measuring system

Measuring method	Light refraction
Measuring Temp min	10 °C

Measuring Temp max	30 °C
Readout	0,2 % 0,001 SG Wort

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	215×94×65 mm
Gross weight	0,35 kg
Net weight	0,30 kg
Shipping method	Parcel service
Packaging component - by weight - paper	888,76 mg
Packaging component - by weight - cardboard	49,9797 g
Packaging component - by weight - plastics	1,41156 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Optional Services

ISO Calibration	961-290
-----------------	---------

Functions

Standard

KERN ORA 3AB

Refractive index measurement for laboratories and the industry

