

KERN OZL 468C825

The flexible, affordable all-rounder with zoom function as a digital solution for schools, training companies, inspection authorities and laboratories.



Article codes, product groups & Status information

Model article-code	OZL 468C825
Model-Serie	OZL-S
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Main article
EAN code (Model)	4045761228719
Planned start date	2019-2HJ
Customs tariff number	90118000
CE Label	Yes
Country of Origin	CN

Technical data - Primary

Optical System	Greenough
Tube	Trinocular
Objectives	0,7x - 4,5x
Illumination	3W LED (transmitted light)

Technical data - Construction

System	Upright
Tube angle	45°
Tube rotation	Yes
Stand type	Mechanical
Stand area	290×240 mm
Holder diameter	76 mm

Eye distance (max)	75 mm
Eye distance (min)	55 mm
Adjustable dioptre	Both-sided
Focusing mechanics	Coarse drive
Visual field min	4,4 mm
Visual field max	28,6 mm
Working distance	105 mm
Sample height max	140 mm
Max resolution	3,7 mm
Dimensions (W×D×H)	300×240×490 mm

Technical data - Power supply

Power supply kind	Built-in power supply
Input Voltage	100 - 240 V
Power Plug	Cold device cable
Protection Class	Protection class I - protective conductor
Cable attachment	Removable
Interfaces	USB 2.0

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	590×440×390 mm
---------------------------------	----------------

KERN OZL 468C825



The flexible, affordable all-rounder with zoom function as a digital solution for schools, training companies, inspection authorities and laboratories.

Gross weight	7 kg
Net weight	5,0 kg
Shipping method	Parcel service
Delivery time	1 d

Technical data - Verification

WEEE	Yes
------	-----

Technical data - Illumination

Illumination directions	Transmitted light
Transmitted light type	LED
Transmitted light intensity	3W
Incident light type	LED
Incident light intensity	3W
Dimmable	Incident light + Transmitted light, separate

Functions

Standard

