

Scoop AISI 304



The stainless steel scoop can be used for taking samples of powders, granulates and pastes. Made of stainless steel AISI 304 (1.4301), the material is sterilisable and autoclavable. The rustproof stainless steel scoop has a smooth, high-gloss polished surface and is food-safe.

The tubular shape makes it easy to immerse the scoop in the medium being sampled. The high sides prevent the sample from trickling out. The round handle offers a good grip. The sturdy scoop is ideal for chemical, pharmaceutical, animal feed and food applications.

- Stainless steel AISI 304 (1.4301)
- Sterilisable
- Robust

Scoop AISI 304

Capacity ml	Dia. mm	Item no.	
120	50	5370-0120	N3
450	75	5370-0450	N3
1000	100	5370-1000	N3
1500	120	5370-1500	N3
2500	150	5370-2000	N3



FoodScoop



Very stable hand scoop for various applications, particularly suited for the food industry. Seamless connection of scoop and hollow handle, therefore without crevices and edges. Especially smooth, highly polished surface.

- Especially sturdy model
- Stainless steel AISI 304 (1.4301)



FoodScoop

Capacity ml	Total length mm	L x W mm	Item no.	
100	215	120 x 60	5371-0100	N3
200	230	140 x 75	5371-0200	N3
250	245	155 x 85	5371-0250	N3
350	290	180 x 100	5371-0350	N3
500	350	220 x 130	5371-0500	N3

