

## Hose tubing connectors



Conical, straight hose connectors for various hose sizes.

- PP
- Sterilisable

## Reducing hose connectors



Reducing nozzles for the reduction from bigger to small hoses.

- PP
- Sterilisable

## Hose couplings



Attachable, conical hose couplings for pressureless operation.


- PE

### Hose tubing connectors, straight, conical



| For dia. mm | Pieces per pack | Item no.     |
|-------------|-----------------|--------------|
| 3-5         | 10              | 8700-0406 N3 |
| 5-7         | 10              | 8700-0608 N3 |
| 7-10        | 10              | 8700-0810 N3 |
| 9-12        | 10              | 8700-1012 N3 |
| 11-14       | 10              | 8700-1214 N3 |
| 13-15       | 10              | 8700-1416 N3 |

### Reducing hose connectors



| For dia. mm | Pieces per pack | Item no.     |
|-------------|-----------------|--------------|
| 4-8/8-12    | 10              | 8703-0408 N3 |
| 4-8/12-16   | 10              | 8703-0412 N3 |
| 8-12/12-16  | 10              | 8703-0812 N3 |

### Hose couplings



| For dia. mm | Pieces per pack | Item no.     |
|-------------|-----------------|--------------|
| 3-5         | 10              | 8704-0406 N3 |
| 5-7         | 10              | 8704-0608 N3 |
| 7-10        | 10              | 8704-0810 N3 |
| 9-12        | 10              | 8704-1012 N3 |
| 11-14       | 10              | 8704-1214 N3 |
| 13-16       | 10              | 8704-1416 N3 |