


PRODUCT DATA SHEET

Delicate Scissors (Iris Scissors), curved, delicate pattern, sharp/sharp, 110 mm (4 1/4") (BC061R)



Attribute Name	Attribute Value
Product Search Type	Delicate Scissors
Altern. Product Name	Iris Scissors
Tungsten Carbide	no
Geometry	curved
Form	delicate pattern
Length	110.00 mm
Length (Inch)	4 1/4"
Tip profile	sharp/sharp
Discipline	Cardiac / Thoracic Surgery, Dental Surgery, Ear-Nose-Throat Surgery, General / Visceral Surgery, Gynaecology, Hand / Foot Surgery, Vascular Surgery, Urology, Spine Surgery, Plastic Surgery, Pediatric Surgery, Orthopaedics / Traumatology, Neurosurgery, Maxillo-Facial Surgery
Indication	Dental Implantology, Dental Surgery
UDI Direct Marking	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No

Attribute Name	Attribute Value
Quality standard	Classic
Shape working end	curved
Handle design	Ring handle
Product category	Dissecting Scissors