

schuett **petriturn**

Turntables for Petri dishes & Inoculating hooks



schuett **biotec**.de



schuett petriturn-E with foot-switch and turntable-attachment for Petri dishes up to Ø 150 mm (optional)



Inoculating hook made of stainless steel (Drigalski-hook)



Inoculating hook made of glass (Drigalski-hook)

Easy & Uniform

Petri dishes turntables

For uniform plating of liquids on agar-surfaces in Petri dishes. Easy to handle, tireless.

schuett petriturn-E Electrically driven Petri dish turntable

With electrically drive and sensor. Body made of stainless-steel. Handling of the turntable contactless by sensor. For Petri dishes up to Ø 100 mm (up to Ø 150 mm optional).

Continuously adjustable, constant speed of 10-100 rpm and continuously adjustable operating time of 3 s to 2 min or permanent running.

Optional connection of a foot-switch.

Turntables for Petri dishes of other diameters available on request.

Comfortable

schuett petriturn-M Manual operation

Powerful and long-running due to high mass of turntable.

Useable on both sides, that means two work spaces in one device by easiest turn-over. For Petri dishes with Ø 100 and Ø 150 mm diameter.



schuett petriturn-M with inoculating hook (optional)

Technical data

petriturn-E

Dimensions (Ø x h)	160 x 70 mm
Weight	approx. 1 kg
Wall power supply	12 V DC 5 VA
Speed continuously, adjustable	up to 10-100 rpm
Operating time, continuously adjustable	3 s - 2 min / permanent running

petriturn-M

Dimensions (Ø x h)	160 x 45 mm
Weight	approx. 2.3 kg

Ordering Information

	Cat.-No.
<i>schuett petriturn-E, electrically driven, 230 V</i>	3.361 252
<i>schuett petriturn-E, electrically driven, 115 V</i>	3.361 262
<i>schuett petriturn-M, manual operation</i>	3.361 102

Accessories

<i>Foot-switch for schuett petriturn-E</i>	3.361 302
<i>Turntable-attachment for Petri dishes up to Ø 150 mm</i>	3.361 402

Inoculating hook

<i>Inoculating hook made of stainless steel (Drigalski-hook), 180 mm length</i>	3.362 152
<i>Inoculating hook made of glass (Drigalski-hook), 145 mm length</i>	3.362 162

schuett[®]biotec.de

Rudolf-Wissell-Straße 13, D-37079 Göttingen, Germany
Fon +49 (0) 551/50 410-0, Fax +49 (0) 551/50 410-99
E-Mail: info@schuett-biotec.de
www.schuett-biotec.de