


PRODUCT BROCHURE

KERN 328-04

Date: 20.04.2025

ENGLISH



Milligram weights, F1, 10 mg, Platelet, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Milligram weights, F1, 10 mg, Platelet, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	328-04
Model-Serie	328
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761240346
Predecessor Model	328-04-2019a
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 mg
Tolerance (OIML)	0,025 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	7,5×6×3 mm
--------------------	------------

Material stainless steel polished

Technical data - Measuring system

Nominal value	10 mg
OIML class	F1
OIML conform	yes
Tolerance (OIML)	0,025 mg
Density (kg/m ³)	7960
Volume (cm ³)	0,0012

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	7.5×6×3 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkKS Calibration	962-454
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
338-090-200	Wooden Box for Set of Weights KERN 338-090-200
347-009-400	Plastic Weight Box KERN 347-009-400
317-009-600	Aluminium Protective Box KERN 317-009-600

SERVICES

Model	Description
962-454R	DAkkS Recalibration KERN 962-454R
920-30	Marking with stroke number KERN 920-30
952-454	Verification KERN 952-454
962-454	DAkkS Calibration KERN 962-454
921-454	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-454
920-301	Marking with stroke number (Last digits of serial number) KERN 920-301

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
328-01	1 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
328-02	2 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
328-03	5 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
328-04	10 mg	0,025 mg	317-009-600	347-009-400	338-090-200	-
328-05	20 mg	0,03 mg	317-009-600	347-009-400	338-090-200	-
328-06	50 mg	0,04 mg	317-009-600	347-009-400	338-090-200	-
328-07	100 mg	0,05 mg	317-009-600	347-009-400	338-090-200	-
328-08	200 mg	0,06 mg	317-009-600	347-009-400	338-090-200	-
328-09	500 mg	0,08 mg	317-009-600	347-009-400	338-090-200	-

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood	Weight set case
328-226	1 mg - 500 mg	-	313-010-600	-	-	in aluminium case
328-31	1 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
328-32	2 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
328-33	5 mg	0,02 mg	317-009-600	347-009-400	338-090-200	-
328-34	10 mg	0,025 mg	317-009-600	347-009-400	338-090-200	-
328-35	20 mg	0,03 mg	317-009-600	347-009-400	338-090-200	-
328-36	50 mg	0,04 mg	317-009-600	347-009-400	338-090-200	-
328-37	100 mg	0,05 mg	317-009-600	347-009-400	338-090-200	-
328-38	200 mg	0,06 mg	317-009-600	347-009-400	338-090-200	-
328-39	500 mg	0,08 mg	317-009-600	347-009-400	338-090-200	-
328-22	1 mg - 500 mg	-	-	328-222-420	-	in plastic case