Flasks, conical, wide neck

PYREX®

Actions Available



- Complies with DIN 12385
- Erlenmeyer pattern
- Manufactured from Pyrex® borosilicate glass
- Excellent thermal performance
- High resistance to chemical attack
- White enamel graduations and marking spot

Item Code	Capacity ml	Body OD mm	Height mm	Neck ID mm
1140/02M	50	48	85	30
1140/04D	100	64	110	30
1140/08D	250	85	140	46
1140/10D	500	105	175	46
1140/14D	1000	131	220	46
1140/18D	20000	153	275	66