MACHEREY-NAGEL

Data sheet

Filter paper

Parameter*	Value*
Grade	MN 621
Basis weight	130 g/m²
Thickness	0.26 mm / 10.2 mil
Filtration speed	43 s
Surface	Smooth
Color	Brownish
Material	Cellulose
Characteristics / Application	Wet strengthened, medium fast filtration, e.g. for soil analysis

*These data are typical values and must not be regarded as specification.

TÜVRheinland CERTIFIED
Management System
EN ISO 13485:2016
ISO 9001:2015
www.tuv.com
ID 0000056401

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany
DE +49 24 21 969-0 info@mn-net.com
FR +33 388 68 22 68 sales-fr@mn-net.com
CH +41 62 388 55 00 sales-ch@mn-net.com
US +1 888 321 62 24 sales-us@mn-net.com